

## Skladový list PromegaExpress™- PEx East Port Praha, UNIVERZITNÍ KAMPUS, Brno, ČR

QTY	CAT#	PRODUCT	LOCATION
<b>RT &amp; PCR Components</b>			
10	M3171	GoTaq® DNA Polymerase, 100u	A1-A5 & B1-B5
5	M5001	GoTaq® Hot Start Polymerase, 100u	C1-C5
3	M7741	<i>Pfu</i> DNA Polymerase, 100u	D1-D3
2	A6001	GoTaq® qPCR Master Mix, 200 reactions	GROUND
4	C1181	Random Primers, 20µg	E1-E4
3	M3681	M-MLV Reverse Transcriptase, RNase H Minus, Point Mutant, 2 500u	F1-F3
3	A3802	ImProm-II™ Reverse Transcriptase, 100 reactions	F4-F6
4	N2511	Recombinant RNasin® RNase Inhibitor, 2 500u	G1-G4
3	N2515	Recombinant RNasin® RNase Inhibitor, 10 000u	G5-G7
4	N2611	RNasin® Plus RNase Inhibitor, 2 500u	H1-H4
3	N2615	RNasin® Plus RNase Inhibitor, 10 000u	H5-H7
<b>Cloning</b>			
2	A1360	pGEM®-T Easy Vector System I	I1-I2
3	M1801	T4 DNA Ligase, 100u	I3-I5
<b>Protein analysis</b>			
5	V5111	Seq Grade Modified Trypsin, 100µg (5 x 20µg)	J1-J5
<b>Restriction enzymes</b>			
1	R6161	Xho I, 3 000u	K1
1	R6041	Hind III, 5 000u	K2
1	R6021	BamH I, 2 500u	K3
1	R6011	EcoR I, 5 000u	K4
1	R6351	EcoR V, 2 000u	K5
1	R6111	Pst I, 3 000u	L1
1	R6801	Nde I, 500u	L2
1	R6513	Nco I, 200u	L3

Dana Lišková, Ph.D., +420 602 570 013

Ondřej Ptáček, Ph.D., +420 602 200 982

Správce mrazáku: Mgr. Tomáš Kašparovský, Ph.D., +420 728 553 721

East Port Praha s.r.o., [eastport@eastport.cz](mailto:eastport@eastport.cz), +420 235 318 177, 800 100 529, [www.eastport.cz/pex](http://www.eastport.cz/pex)